
PATIENT CONVEYOR KIT





A quick deployment patient conveyor in three sections for decontamination procedures and patient handling that requires minimal personnel to deploy.

The patient transfer board is designed to work with the conveyor to allow efficient and effective treatment of nonambulatory patients. The transfer board design folds for storage and transportation and locks into a sturdy platform when in use.

FEATURES

- Designed specifically to enable rapid responders to effect NBC decontamination processes within minutes.
- Lightweight geodesic framework manufactured from high strength powder coated aircraft alloy.
- New transfer board design folds for storage and transportation and locks into a sturdy platform when in use.
- Includes locking plates to fix the conveyor in the open position and allow multiple conveyors to be connected together.
- Storage and transportation bags included as standard.
- Three section design for use with the MFC Mass Decontamination Shower, allowing the central section of the shower to be closed when in use

TECHNICAL DATA

Patient Conveyor & Transfer Board

10

Product Code	DSA015 6.8m Conveyor (total) & Transfer Board	
Product Name		
Conveyor Dimensions	<u>Open</u>	Closed
- Length (cm)	220 , 240, 220 (Section 1, 2 & 3)	96, 96, 96 (Section 1, 2 & 3)
- Width (cm)	50	15
- Height (cm)	76	79
Conveyor Weight (kg)	Total 90 (each section approx 30kg)	
Transfer Board Dimensions	<u>Open</u>	Closed
- Length (cm)	211	106
- Width (cm)	59	59
- Height (cm)	3	7

3: www.mfc-international.com

Board Weight (kg)

!: +44 (0) 1443 433 075

MFC International